

SADS UK Cabinet order form

Name of organisation ordering this cabinet:

Delivery address & Tel #

Invoice address:

.....

Your Reference/purchase order :

Signature Printed name Date.....

Tick if VAT exempt Certificate enclosed Ambulance Service contact:.....

Payment:-30 days from invoice. Account name SADS UK Sort code 60-03-25 Account# 77056965

Stainless steel Cabinets	£	Qty
Key coded locking Cabinet, colour yellow (RAL 1016)	350.00	
Non-locking Cabinet colour yellow (RAL 1016)	350.00	
Polycarbonate Cabinets 4000 Series		
Non-Locking Cabinet colour yellow – Ex display cabinet	350.00	
Key coded Locking Cabinet Colour yellow – Ex display cabinet	350.00	
Tough Plastic GRP		
Non-Locking Cabinet partially painted	330.00	
Key coded locking Cabinet partially painted	330.00	
And		
Non-Locking Fully painted	350.00	
Key coded locking Fully painted	350.00	

The above prices do not include VAT or P&P

All cabinets are provided with **battery** powered alarm.

All heated cabinet requires a 240v supply. Installation must be undertaken by a suitably qualified person, this cost is not included in the above. Our cabinets use a double insulated 240v without earth wire. Power consumption approx. 50w

Orders to be sent:

SADS UK, Suite 1a Churchill House, West Horndon, Essex CM13 3XD email: technical@sadsuk.org Tel: 01277 811215

Important information

Always dial 999 in an emergency. The defibrillator can be used to save a live in a cardiac arrest emergency whilst the ambulance service is on its way.

Keycode locks. It is important to advise your ambulance service of the new defibrillator site. It must be registered on the BHF Circuit, SADS UK will help you with this. When a 999 call is made the defibrillator can then be located and the keycode provided to allow access to the defibrillator.

Put the guardians contact details on the back of your defibrillator. If the defibrillator is taken to an emergency it may go to hospital with the patient, having your details on the back will help it to be returned to you. Add the cabinet keycode lock number to these details so that the defibrillator can be put back in the cabinet after an incident.

Visual checks of the defib status indicator must be carried out regularly and consumables replaced when necessary. **Please refer to SADS UK checklist. Regular opening of the cabinet and checking the defibrillator status is essential.**

Cabinet postcode - Ensure you include the cabinet postcode using indelible ink on the outside of your cabinet so that the rescuer can provide this location to the 999 call handler. If the ambulance service wants to use an additional reference they will advise you.